

HOW TO USE

TURN ON/OFF

TURN ON AND OFF WAND

Press and hold the power button for 3 seconds until the display lights up. Do the same to TURN OFF.

TURN ON HEAT

Press the temperature (°C) button to TURN ON the HEAT.



TEMPERATURE CONTROL

INCREASE TEMPERATURE

The temperature will start off at 36 degrees, to increase the temperature click the °C button, each click will increase temperature by 1 degree. Maximum temperature 53 degrees.

DECREASE TEMPERATURE

Press and hold the °C button until it reaches 53 degrees and one more click will restart the temperature to 36 degrees.

PROTECTION & DISPLAY

TEMPERATURE DISPLAY

Each individual finger will read a temperature on the top of the display, this will vary +-3 degrees depending on the environmental conditions.

FINGER TIPS

Only the tips of the metal fingers, past the white band, will get hot.

PROTECTIVE COVER & CASE

Used to protect the metal fingers when not in use. Please remove when in use. Store wand in case provided like in the picture.



EXTRA FEATURES

MICROVIBRATION

Press the button on the top right to turn on, starts on the lowest setting of the 3 settings available, simply click through each one. Click after the highest setting will turn the vibration off.

INFRARED LIGHT

Press the button on the top left to turn on and press again to turn off.

PAUSE

Vibration and light therapy simply click the power button once. This won't pause the heat.

CHARGING

At the bottom of the device is the USB port. Use the USB charging cable supplied and leave on CHARGE for 3-4 HOURS until the lights stop flashing red and green.

When FULLY CHARGED the light will stay GREEN.

Battery lasts for 2-3 hours depending on the temperature and other settings used.

WHERE TO USE

Gently Massage the heated metal fingers into the tight muscles on the body, of both horses and people, to relax them. It is also suitable for other animals such as dogs.

Not just for tight muscles, the wand helps aid relaxation and reduce stress.

HORSES

PEOPLE



NECK



BACK OF HIND LEG



AROUND SHOULDER
BLADE



LOWER BACK



POLL



SHOULDER



BACK OF ARM



NECK AND
TOP OF SHOULDER



BACK



QUARTERS

EXAMPLES OF THE
BEST PLACES
TO USE
THE THERMAL WAND...

TECHNICAL SPECIFICATIONS

Charging Voltage:	DC5V 1A
Charging Time:	3-4 Hours
Use Power:	6.5W
Standby Time:	Around 40 Days
Vibration Frequency:	10000rMb+ 10%
Battery Volume:	1600mA/H:
Working Time:	2.30-3.30 Hours
Product Weight:	380g
Product Size:	20cm*6.9cm*3.4cm
Product Safety Certificates:	FCC ROHS CE: